S E R V I C E N O T E

SUPERSEDES: NONE

4156B Semiconductor Parameter Analyzer 04155-66541/04155-69541 CPU Assembly

Serial Numbers: See below.

CPU board replacement to prevent startup failures

Duplicate Service Notes:

4155B-05 Serial Numbers: JP10E00240/JP10E00279 4156B-05 Serial Numbers: JP10E00290/JP10E00331

To be Performed by: Agilent-qualified personnel

Parts Required: 04155-69541

Situation:

The unit does not start up and the following message may be displayed:

INSERT UPDATE MASTER DISK, THEN CYCLE POWER.

Continued

DATE: January 1999

ADMINISTRATIVE INFORMATION

SERVICE NOTE CLASSIFICATION:		
MODIFICATION RECOMMENDED		
ACTION CATEGORY:	☐ IMMEDIATELY ☐ ON SPECIFIED FAILURE ■ AGREEABLE TIME	STANDARDS: LABOR 0.5 Hours
LOCATION CATEGORY:	☐ CUSTOMER INSTALLABLE☐ ON-SITE☐ SERVICE CENTER	SERVICE RETURN USED RETURN PARTS: SCRAP SEE TEXT SEE TEXT SEE TEXT
AVAILABILITY:	PRODUCT'S SUPPORT LIFE	AGILENT RESPONSIBLE UNTIL: January 2001
AUTHOR: MO	ENTITY: 3300	ADDITIONAL INFORMATION:

© 1998 AGILENT TECHNOLOGIES PRINTED IN U.S.A.



Page 2 Service Note 4156B-05

Solution/Action:

Replace the CPU board assembly (P/N 04155-69541). Do NOT try to reload the firmware.

The above problem is caused by defective flash ROMs on the CPU board. Therefore, the problem may recur if you reload the firmware.

Flash ROMs that may cause the problem are identified by the revision and date code information on the ROMs. Remove the shield plates of the CPU board and check the information.

AM29F040 ... Good

AM29F040B ... Check the date code ... 9752GBM J (9752 is a date code of the defective lots.) The date codes 9804 and later are good.

When an 4155B or 4156B in the above serial number bands is returned to an Agilent service center, check the flash ROM revision and date code. If the flash ROMs are in a defective lot, replace it with a good board (AM29F040, or AM29F040B with date code 9804 or later). If you find a replacement CPU board in a defective lot, do not use it and return it.

Perform this service note in "Agreeable Time" category for the above serial number bands even if the symptom is not observed.