# M9018A-04

# S E R V I C E N O T E

Supersedes: NONE

## M9018A 18 slot PXIe Chassis

**Serial Numbers: TW51060101-TW51510103** 

## M9018A Serial Number Does Not Display in Agilent Connection Expert

**Parts Required:** 

P/N Description Qty.

Part Number	Description	Quantity
M9018-69103	18-slot PXIe Refurbished Exchange Unit Assembly (exclude fan	1
	assembly)	

### ADMINISTRATIVE INFORMATION

SERVICE NOTE CLASSIFICATION:					
MODIFICATION RECOMMENDED					
ACTION √ ON SPECIFIED FAILURE CATEGORY: [[]] AGREEABLE TIME	STANDARDS  LABOR: 1.0 Hours				
LOCATION [[]] CUSTOMER INSTALLABLE CATEGORY: [[]] ON-SITE (active On-site coverage only)  √ SERVICE CENTER [[]] CHANNEL PARTNER	SERVICE √RETURN INVENTORY: [[]] SCRAP [[]] SEE TEXT USED √ RETURN PARTS: [[]] SCRAP [[]] SEE TEXT				
AVAILABILITY: PRODUCT'S SUPPORT LIFE	NO CHARGE AVAILABLE UNTIL: Always				
[[]] Calibration Required √ Calibration NOT Required	PRODUCT LINE: BL AUTHOR: TP				
ADDITIONAL INFORMATION:					

© AGILENT TECHNOLOGIES, INC. 2013 PRINTED IN U.S.A.



February 13, 2013

Rev. 21


Page 2 of 2 M9018A-04

#### **Situation:**

It has been found for the affected M9018A described above with the configuration setting of:

- 1. "2-Link Configuration: 2x8" and using the Agilent M9036A PXIe Embedded Controller
- 2. Or "4-Link Configuration: 4x4" and using 3<sup>rd</sup> party embedded controller

That the M9018A chassis serial number does not display when opening Agilent Connection Expert. See picture below:



#### **Solution/Action:**

If the symptom is observed, contact Agilent Service for repair.

### **Revision History:**

	Revision Number	Date	Author	Reason For Change
1.0	0	2/7/2013	TP	As published